



# Total Lift Bed™

Early Mobility at the Push of a Button™

When healthy, we rely on regular movement to promote blood flow, muscle growth, and a myriad of other physical benefits. Unfortunately, not everyone is able to move freely on their own. In critically ill patients, prolonged bed rest and immobility can lead to severe physical deconditioning, significant loss of muscle mass, and increased muscle weakness that can result in extended time spent in both the ICU and the hospital. Early mobilization can help alleviate the adverse, long-term physical and psychological consequences of prolonged immobility—and automated mobilization solutions can reduce risks associated with manual lifting to improve outcomes for patients and caregivers alike.

### **Early Mobility at the push of a button™**

Using only one or two caregivers, the Total Lift Bed®\* can progressively move patients to a full standing and weight bearing position in order to aid rehabilitation early in the acute phase of critical illness.

### **Compliant with early mobility protocols**

Without leaving the bed, patients who need a graded transition from bed rest to supporting their full body weight can be progressively moved to the standing position as early and frequently as tolerated.

### **Patented foot Lifter footboard**

When the tilt function is activated, the patented Foot Lifter™ moves towards the patient's feet. When contact is made, the bed begins gently tilting to the desired angle or weight bearing ratio, while maintaining the patient's position to help avoid the need for repositioning. The Foot Lifter will automatically adjust to become even with the floor for patient egress from the bed if appropriate.

### **Safety during movement and while standing**

A flexible strap system helps secure the patient during movement to the upright position.



**C-ARM Compatible**

Patients from all around the world are benefiting from the Total Lift Bed



## Chair Mode

The Chair mode on the Total Lift Bed Is Unique in itself. During the the movement the bed uses it's distinctive design to dip the bed back at the hips to help prevent sliding and provide a more comfortable seating position with less pressure on the sacrum.

## Trendelenburg

The newest version of the Total Lift Bed now does Trendelenburg posiiton too.

## Total Drive

Total Lift Bed now comes with the option to add our Patented new power drive system. Unique in that you don't have to raise and lower the system seperatly because the patented design allows sideways and front and back movement even without power. Speed controlled and easy to use.

## Mattress Options

Total Lift Bed now comes with more Mattress Options. Even the standard comes with a stage 4 alternating pressure low air loss mattress. Nwo you can also order turn assist, rotational, pulsation and self adjsuting mattresses.





3520 NW 56th St.  
Fort Lauderdale, FL 33309

tel. 877-678-6781

## Studies



Piedmont Makes Protocol for Ambulating Patients with Aortic Ballon Pump.



University Of Texas shows a substantial decrease in length of Stay using Total Lift Bed.

